## Erkoplast PLA-R rosa

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### EN: Instructions for use

**Intended use:** Thermoforming sheets and accessories for the fabrication of dental and orthodontic appliances for the treatment of all patient target groups.

**Field of application:** Erkoplast PLA-R is exclusively suitable for the fabrication of base plates and bite plates. Fabrication, see thermoforming technique brochure and videos, [www.erkodent.com](http://www.erkodent.com)  

**Contraindications:** In case of already known allergic reaction to plastics a specific allergy to Erkoplast PLA-R (PLA) should be excluded.

**Precautionary measures:** Erkoplast PLA-R is only recommended for the above mentioned devices. Observe storage instructions.

**Attention:** The use is limited to the production of dental and orthodontic appliances. Use is subject to the responsibility of the therapist. Use only on the recommendation or prescription of an authorised therapist. For single use only (thermoforming). Allergic reactions are unlikely, but cannot be excluded. Incorrect processing can lead to cracks, fractures and sharp edges. There is a risk of aspiration in the case of small breakages. Serious incidents in connection with Erkoplast PLA-R must be reported by the user to the competent authority and to the manufacturer.

Erkoplast PLA-R (approx. values): Density 1.24 g/cm<sup>3</sup>, modulus of elasticity 3500 MPa, elongation at break 5 %, water absorption 0.1 %, temp. resistance 55 °C

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