

Product data sheet Aton-OF 80



1. Manufacturer

Trade name:	Aton-OF 80
Purpose:	Production of impressions
Manufacturer:	Erkodent Erich Kopp GmbH Siemensstraße 3 72285 Pfalzgrafenweiler Germany Phone: +49 7445 8501-0

2. Application

Kneadable addition-curing silicone for precise impressions of all kinds.

3. Properties

Technical features:

Properties	Test method	Value
Form	ISO 4823:2015	kneadable
Color	-	White/yellow
Smell	-	characteristic
Density	ISO 1183	1.65 g/cm ³
Hardness according to Shore A	ISO 868	± 80
Best before	-	36 months

4. General information

Storage Notes:

Store away from sunlight and in a dry place.

Recommended storage temperature: 5°C – 35°C

Recommended processing temperature: 15 – 25°C

Cleaning and care instructions:

Do not come into contact with latex gloves, disinfectants, foot creams or similar.

To the best of our knowledge, the information on this product data sheet corresponds to our findings at the time of going to press. The information is not a guaranteed assurance of product characteristics and does not constitute a contractual understanding of the law. Input errors and errors excepted.