

Product datasheet Atran-OF



1. Manufacturer

Trade name:	Atran-OF
Purpose:	Production of impressions
Manufacturer:	Erkodent Erich Kopp GmbH Siemensstraße 3 72285 Pfalzgrafenweiler Germany Phone: +49 7445 8501-0

2. Application

Paste-like addition-curing silicone for the production of detailed impressions on foot models.

3. Properties

Technical features:

Properties	Test method	Value
Form	ISO 4823:2015	pasty, slightly flowing
Color	-	White/Brown
Smell	-	characteristic
Density	ISO 1183	1.23 g/cm ³
Hardness according to Shore A	ISO 868	± 41
Best before	-	24 months

4. General information

Storage Notes:

Store away from sunlight and in a dry place.
Recommended storage temperature: 5°C – 35°C
Recommended processing temperature: 15 – 25°C

Cleaning and care instructions:

Do not bring into contact with latex, disinfectants, foot creams or similar.

To the best of our knowledge, the information on this product data sheet corresponds to our findings at the time of going to press. The information is not a guaranteed assurance of product characteristics and does not constitute a contractual understanding of the law. Input errors and errors excepted.