

Product datasheet Erkosit liquid



1. Manufacturer information

Trade name:	Erkosit liquid
Intended use:	Fabrication of temporary appliances
Manufacturer:	Erkodent Erich Kopp GmbH Siemensstraße 3 72285 Pfalzgrafenweiler Germany Tel.: +49 7445 8501-0

2. Intended use

Self-polymerizing resin for the fabrication of temporary appliances. Especially suited for the durable combination with Erkodur using PV-Primer.

3. Properties

Technical properties:

Properties	Guideline	Value
Form	ISO 4823:2015	liquid
Colour	-	clear
Odour	-	characteristic
Density	ISO 1183	0,97 g/cm ³
Flashpoint	-	10°C
Solvent content	-	100 %
Minimum durability	-	36 months

4. General information

Storage instructions:

Keep away from sunlight. Keep dry.

Recommended storage temperature: 5°C – 35°C

Recommended working temperature: 15°C – 25°C

The information on this product data sheet corresponds to the best of our knowledge at the time of printing. The information is not guaranteed assurance of product properties and does not constitute a contractual legal understanding. Input errors and mistakes reserved.